



News – For Immediate Release

Date: July 29, 2015

CPAC Announces New COX RapidHeat Sterilizer Design

The COX Rapid Heat Model 6000 - the industry leader in high-velocity hot air sterilization featuring FDA approved 6 & 12 minute cycles introduces its new look in functional cabinet design.

This HVHA (High-Velocity High Temperature) sterilization system does not require the use of any water, steam or toxic chemicals. Drying cycles are eliminated in HVHA systems and with the absence of water, steam and chemicals there is no corrosion, pitting or dulling of instruments.

The COX Rapid Heat Sterilizer has been recently proven by an independent 3rd party to average 84% less energy per sterilization cycle than competitive steam sterilizers, and can sterilize up to 150 instruments per hour in pouched cassettes. The fast turnaround of the COX Rapid Heat Sterilizer eliminates the temptation of medical staff to bypass or shorten cycles that can jeopardize instrument sterility.

In addition to rapid dry heat sterilizers, the company offers a full complement of instrument sterilization pouches, instrument organizers, cassettes, holding racks, spore test kits, incubators and indicator strips.

For more information contact: Diane Baker, dbaker@cpac.com - Toll Free 800-828-6011, extension 1372.
