



News – For January 2019 Release

Date: 01-07-2019

High-Capacity RapidHeat Tabletop Sterilizer Now Available

After a period of lengthy development, CPAC Equipment, Inc. is pleased to announce the February 2019 rollout of their new RapidHeat™ Pro11 High-Capacity Tabletop Sterilizer.

With over 700% (7 times) more capacity than its COX RapidHeat™ predecessor, the Pro11 has the industry's fastest "complete" cycle turnaround. Unlike steam sterilization there is no FDA required drying cycle that can be short-cut to speed up the process. The Pro11's complete full cycle from start to finish including warm-up is 12 minutes Unwrapped or 18 minutes Wrapped. The Time/Temperature profile established for each cycle has been documented to exceed a 12 Log₁₀ bacterial spore kill. This documentation is based on a sterilizer capacity of four (4) instrument trays, each containing a load of over 1.75 lbs. representing a total sterilizer capacity of approximately 140 medical or dental instruments.

And since no water or steam is used, independent tests have confirmed that there is less maintenance expense with each Pro11 cycle consuming 85% less energy than comparable steam sterilizers. These tests have also confirmed that the RapidHeat™ waterless sterilization environment does not contribute to the instrument corrosion and shortened instrument life that has been proven to be characteristic of steam sterilization.

For more information: Call Toll Free 800-828-6011, ext. 1372 or E-mail dbaker@cpac.com.